

# Sanapur®

Recommended by physiotherapists.

## Slatted frame 2M - flat motor frame with 2 motors

### For an optimal lying position

The SANAPUR slatted frame completes the SANAPUR sleeping system so that the mattress can optimally develop its properties. By means of the height-adjustable slats, the slatted frame is individually adapted to your ergonomic requirements.



*The slatted frames can be used in almost all common bed surrounds.*



### Individual comfort

The **height-adjustable ledges** (anthracite) in the shoulder, lordosis and pelvic zones are 5-fold adjustable. Starting from the central 0-position, the slats can be lowered up to two positions (more relief) and/or raised (more support).

### Product description

#### Solid beech outer frame

- 70 mm high / 20 mm thick, stable and distortion-free

#### Spring laminated beech wood slats (6-fold glued)

- 12 fixed slats in beech nature
- 16 height-adjustable anthracite slats (shoulder, lordosis and pelvic zones)

#### caps

- 28 protruding caps made of high-quality and tested plastic (REACH)

#### dimensions

- Available in: Lengths 190/200/210/220 x widths 80/90/100/120/140 cm (nominal dimensions)
- Frame height approx. 95 mm, height adjustable slats max. 105 mm
- Real frame dimensions: length 40 mm shorter / width 15 mm narrower than nominal dimension

### Functional description

#### Flat motor integrated in foot area

- Motor flush with the lower edge of the frame, perfect for drawer runners / bed drawers
- Rope pull technology on both sides ensures uniform drive

#### Useful Information

- Included in delivery: cable remote control and 1 mattress holder foot end
- When the upper body is lifted, the headrest is pulled back slightly delayed
- Two additional metal struts stabilize the foot section
- Mains disconnection and emergency lowering
- Frame TÜV certified up to 150 kg load (incl. mattress)

#### Accessories (optional)

- Foot Set (easy assembly)
- Mattress holder (lateral)